

Electrical Impedance Tomography: Image Algorithms and Applications

Andy Adler

Assistant Professor,
SITE, U.Ottawa


Outline

- Electrical Impedance Tomography
- Physics
- Image Reconstruction
- Applications
- Future Work

Electrical Impedance Tomography

- Relatively new medical imaging technique (early 1990's)
- Body Surface Electrodes apply current patterns and measure the resulting voltages
- Distribution of conductivity is calculated

EIT: Block Diagram



EIT: Applications


- EIT can image physiological processes involving movement of conductive fluids and gasses
- Lungs
- Heart / perfusion
- GI tract
- Brain
- Breast

EIT: Advantages

EIT is a relatively low resolution imaging modality, *but*


- Non-invasive
- Non-cumbersome
- Suitable for monitoring
- Underlying technology is low cost

Application: Breathing




Chest images of tidal breathing in normal

Application: Heart Beat



EIT signal in ROI around heart and ECG

Non-invasive



Thresholds for cutaneous perception of electric current vs. frequency and EIT system

EIT: Physics

- Within medium Ω there is \mathbf{E} and \mathbf{J} .

- $$\mathbf{J}_c = \sigma \mathbf{E}$$

- $$\mathbf{J}_d = \epsilon \epsilon_0 \frac{d\mathbf{E}}{dt}$$

$$\mathbf{J} = (\sigma - j\omega\epsilon\epsilon_0)\mathbf{E}$$

EIT: Physics

In the absence of magnetic fields

$$\mathbf{E} = -\nabla V$$

No charge build up in conductive medium

$$\nabla \cdot \mathbf{J} = -\frac{\partial \rho}{\partial t} = 0$$

We have

$$\nabla \cdot (\sigma - j\omega\epsilon\epsilon_0)\nabla V = 0 \quad \text{in } \Omega$$


Image Reconstruction: Static Imaging

Static imaging reconstructs the absolute conductivity from measurements.

Algorithms:

- Iterative (Newton-Raphson)
- Layer Stripping

Block Diagram of Iterative Algorithm



Static Imaging Difficulties


- Extremely sensitive to uncertainties in electrode position
- Ill-conditioned problem
- Numerical instability

Dynamic Imaging

- Calculate change in conductivity distribution from change in measurements
- Inverse problem *linearized*
- Much reduced sensitivity to electrode and hardware errors.
- Very suitable for physiological imaging: lung, heart, GI

Dynamic imaging techniques

- *Backprojection*: Voltages are “projected” along the equipotential lines.



Inverse Techniques

- We can pose dynamic imaging as linear inverse, using a *sensitivity matrix*

$$\mathbf{z}_j = \frac{\mathbf{z}(\sigma_h) - \mathbf{z}(\sigma_h + \delta_j)}{\delta_j}$$

$$\mathbf{z} = \mathbf{H}\Delta\sigma$$

Inverse Techniques

- Classic least-squares inverse

$$\mathbf{z} = \mathbf{H}\mathbf{x}$$

$$\hat{\mathbf{x}} = \left(\mathbf{H}^t \mathbf{H}\right)^{-1} \mathbf{H}^t \mathbf{z}$$

Matrix Techniques

However, problem is:

- ill-conditioned: measurements depend much more on data near electrodes than in centre
- ill-formed: more unknowns than measurements

Regularized Imaging

Given Linear Model:

$$\mathbf{z} = \mathbf{H}\mathbf{x} + \mathbf{n}$$

Maximum A Posteriori (MAP) estimate is:

$$\hat{\mathbf{x}} = \left(\mathbf{H}^t \mathbf{R}_n^{-1} \mathbf{H} + \mu \mathbf{R}_x^{-1} \right)^{-1} \left(\mathbf{H}^t \mathbf{R}_n^{-1} \mathbf{z} + \mu \mathbf{R}_x^{-1} \mathbf{x}_\infty \right)$$

Regularized Imaging

- Parameters \mathbf{W} , \mathbf{Q} , \mathbf{x}_∞ , represent *a priori* statistical knowledge of problem

$$\mathbf{x}_\infty = E[\mathbf{x}]$$

$$\mathbf{R}_x = E[(\mathbf{x} - \mathbf{x}_\infty)^t (\mathbf{x} - \mathbf{x}_\infty)] = E[\mathbf{x}^t \mathbf{x}] - \mathbf{x}_\infty^t \mathbf{x}_\infty$$

$$\mathbf{R}_n = E[\mathbf{n}^t \mathbf{n}] = \begin{bmatrix} \sigma_1^2 & 0 & \dots \\ 0 & \sigma_2^2 & \\ \vdots & & \ddots \end{bmatrix}$$

Advantages of Regularization





- Stabilizes ill-conditioned inverse
- Introduction of *a priori* information
- Control of *resolution-noise* performance trade-off
- MAP inverse justifies the formulation in terms of Bayesian statistics

Regularized Inverse

Parameters:

- \mathbf{W} : models measurement noise
- \mathbf{Q} : penalizes image features which are greater than data supports
- \mathbf{x}_∞ : represents the background conductivity distribution (heart, lungs, etc)
- μ : “hyper-parameter” amount of regularization


Noise – Resolution Tradeoff



F: *Meas:* No Noise *Reconst:* NF= 2.0

G: *Meas:* -3dB SNR *Reconst:* NF= 2.0


Noise – Resolution Tradeoff



F: *Meas:* No Noise *Reconst:* NF= 0.4


G: *Meas:* -3dB SNR *Reconst:* NF= 0.4

Electrode Movement



- Electrodes on chest move during breathing
- Figure shows rough movement pattern of rib cage

Electrode Movement



- FEM of electrical and mechanical properties of thorax to simulate
- Signal from expansion shown to contribute 10-20% of conductivity signal

Applications: Lung Function and Disease


Monitoring of lung function is important in ICU. Issues are:

- Distribution of ventilation
- Ventilation – perfusion match
- Lung Edema
- Flow limitation (obstructive lung disease)

Lung function tests


- Patients with obstructive lung disease (emphysema, bronchitis) tend to dynamically hyperinflate: lower inspiration than expiration.
- Simulated with PEEP (positive end expiratory pressure)

Lung Function: Protocol




- Initially, inflations to 3 kPa
- dog expires to the applied PEEP
- ventilation stopped, and dog expires against PEEP for 15 s

Lung Function



A: image of section thorax due to ΔVL 800 ml

B: image due to a ΔVL 400 where left main stem bronchus was plugged.




- Volume estimates by EIT, Pao and Pes, after step volume increases of 100ml, 500ml, 900ml.
- Note that EIT signal does not display overshoot

Pulmonary Edema

- Pulmonary Edema (lungs filled with fluid) plagues ICU patients for heart disease, accident victims
- Typical monitoring techniques are invasive (Swan-Ganz catheter)
- We looked at EIT ability to monitor level of Edema.

Pulmonary Edema: Results



- Change in lung liquid volume by EIT vs liquid volume instilled

Imaging

- Algorithms to use data from internal electrodes
- How to manage electrode artefacts
- Arbitrary electrode placements
- Image interpretation and error bounds given artefacts and noise

Applications of interest

- Longer term monitoring in ICU
- Small animal lung function
- Contrast agents
- Inhomogeneous ventilation on smaller scale